Dickson WT630 WiZARD2 Ethernet Vaccine Temperature Data Logger PDF



Call our friendly team on +44 (0)1243 558270 Tempcon Instrumentation Ford Lane Business Park Ford West Sussex BN18 0UZ, UK www.tempcon.co.uk



Dickson WT630 WiZARD2 Ethernet Vaccine Temperature Data Logger

Product Images



Short Description

Place these data loggers throughout your facility and receive all data automatically at your PC. No retrieving data from loggers, no hard wiring loggers back to your PC. This is the easiest and most cost-effective way to create a scalable monitoring network.

Description

Place these data loggers throughout your facility and receive all data automatically at your PC. No retrieving data from loggers, no hard wiring loggers back to your PC. This is the easiest and most cost-effective way to create a scalable monitoring network.

- One thermistor probe with attached glycol bottle
- Audio/Visual Alarms
- Email and Text alerts available
- Eliminates the need to change charts and pens
- Vaccine data is wirelessly downloaded to your computer
- Loggers continue to log in the event of a power outage

Additional Information

Brand	Dickson
Channel(s)	1
Typical applications	Life Sciences, Pharmaceutical, Vaccine Storage
Measurements	Temperature